

FOR THE ACTIVE WHEELCHAIR USER



Excellent lateral and frontal stability
Perfect fit for the user who transfers frequently
Ischial offloading and protection

The Vicair Active offers the unique combination of ergonomically shaped foam and Vicair SmartCells. The three rear SmartCell compartments provide a high level of comfort and skin protection, and offload the sitting bones and tailbone area. The foam front compartment makes transfers easy. It also gives a great sense of stability and better wheelchair handling. The cushion is made of comfortably breathable textile, that is washable at 60°C. The Vicair Active has an easy accessible storage pouch for precious belongings (keys, phone, wallet). In short: everything you need to maintain your active lifestyle!

PRODUCT PERFORMANCE

Skin protection / Anti decubitus	▲▲▲▲▲
Pressure redistribution	▲▲▲▲▲
Shear reduction	▲▲▲▲▲
Micro climate control	▲▲▲△△
Impact & vibration damping	▲▲▲▲▲

Positioning

Positioning	▲▲▲▲▲
Lateral stability	▲▲▲▲▲
Frontal stability	▲▲▲▲▲
Sliding prevention	▲▲▲▲▲
Allows manual adjustments	▲▲▲▲▲
Allows easy transfers	▲▲▲▲▲

Maintenance & Safety

Maintenance & Safety	▲▲▲▲△
Durability	▲▲▲▲▲
Reliability	▲▲▲▲▲
Functionality when misused	▲▲▲△△
Low maintenance level	▲▲▲▲△
Level of hygiene	▲▲▲▲△

Functionality & Application

Functionality & Application	▲▲▲▲▲
Comfort	▲▲▲▲▲
Low cushion weight	▲▲▲▲△
Available size range	▲▲▲△△
Maximum user weight	▲▲▲▲▲
Low seat to floor height	▲▲▲▲▲
Auto adjustments	▲▲▲▲▲

USED WITH

Spinal cord injuries, active users with a high transfer frequency

TECHNICAL SPECIFICATIONS

Height	6 and 9 cm
Average weight	700 grams
Max. user weight	200 kg
Machine washable	60°C (when emptied)

CUSHION SIZES

- Standard
- On request

(Sizes >53 cm height 9 only)

