### VICAIR® ADJUSTER 02

### THE OPTIMAL CUSHION FOR WHEELCHAIR USERS WITH BODY ASYMMETRY

#### Follows the body's contours

for equal pressure distribution | support ITcoccyx area | managing pelvic asymmetry



## Individual solution for asymmetry

#### Minimizes sliding and shear forces

# Optimal adjustment to the body's contours for equal pressure distribution

The Adjuster O2 is extremely suitable for use with asymmetric body contours and amputations. It automatically manages small obliquities and is individually adjustable to support more extreme asymmetry. When you're looking for a high comfort level, optimal skin protection and stability: the Adjuster O2 is always a safe choice. The design allows total immersion to equally redistribute pressure and limits sliding forward. 100% machine washable and breathable, for optimal hygiene and perfect temperature and moisture regulation. As a result the Adjuster O2 offers you optimal protection against

PRODUCT PERFORMANCE	
Skin protection / Anti decubitus	***
Pressure redistribution	***
Shear force reduction	***
Microclimate control	***
Impact & vibration damping	***
Positioning	***
Lateral stability	
Frontal stability	***
Sliding reduction	***
Allows manual adjustments	***
Allows easy transfers	
Maintanana C Cafatu	
Maintenance & Safety	•

Maintenance & Safety	
Durability	4444
Reliability	<b>AAAA</b>
Functionality when misused	
Low maintenance level	
Level of hygiene	4444

***
<b>AAAA</b>
<b>AAAA</b>
<b>AAAA</b>
4444
<b>AAAA</b>

#### USED WITH

**CUSHION SIZES** 

Hemiplegia (CVA), amputations, foot propulsion, kyphosis, pelvic obliquity/ tilt/rotation, Spina Bifida, spasms, Fibrodysplasia Ossificans Progressiva (FOP)

TECHNICAL SPECIFICATIONS								
Height	6 and 10 cm							
Average weight	725 grams							
Max. user weight	250 kg							
Machine washable	60°C max.							

	Sta	ndard (Sizes >53 cm height 10 onl								only)		
	On	req	ues	st								
												width
Inch		14	15	16	17	18	19	20	21	22	23	24
	СМ	35	38	40	43	45	48	50	53	55	58	60
14	35											
15	38											
16	40											
17	43											
18	45											
19	48											
20	50											
21	53											
22	55											
23	58											
24	60											