VICAIR® MULTIFUNCTIONAL 02 BACK

THE BEST BACK SUPPORT FOR USERS OF MANUAL TILT-IN-SPACE WHEELCHAIRS AND POWERCHAIRS



support the shoulders | compartment shape provides freedom of movement

Top compartments

Multidimensional lumbar support





PRESSURE ULCER TREATMENT

Supporting different postures

Follows the human spine alignment

support

Optimal pressure redistribution

Excellent pressure redistribution, comfort and support make this back support cushion perfect for use in a large variety of situations. It can be used in several types of back systems and in tilt-inspace chairs with adjustable seat depth. The cushion is easily fixed onto the solid back support by use of hook and loop strips. With the Vicair Back Connector the Multifunctional O2 Back can also be used in tension adjustable backs such as the Vicair Strap Back. Although it automatically adapts to the body's contours, the cushion can also be adjusted for optimal positioning, while the user is still leaning against the cushion. 100% machine washable and breathable, for optimal hygiene and perfect temperaa and maistura regulation

PRODUCT PERFORMANCE	
Skin protection / Anti decubitus	
Pressure redistribution	***
Shear force reduction	
Microclimate control	***
Impact & vibration damping	***
Positioning	
Lateral stability	
Frontal stability	
Sliding reduction	***
Allows manual adjustments	***
Allows easy transfers	

Maintenance & Safety	***
Durability	4444
Reliability	AAAA
Functionality when misused	
Low maintenance level	***
Level of hygiene	AAAA

Functionality & Application	
Comfort	AAAA
Low cushion weight	AAAA
Available size range	AAAA
Maximum user weight	AAAA
Low seat to floor height	n.a.
Auto adjustments	

USED WITH

Neuromuscular diseases, Scoliosis, Kyphosis, Lordosis, Gibbus, Quadriplegia, pain relief, in a geriatric environment

TECHNICAL SPECIFICATIONS		
Length	36-60 cm	
Average weight	975 grams	
Max. user weight	150 kg	
Machine washable	60°C max.	
Loss of seat depth	3 cm	

CUSHION SIZES	
Standard	
On request	

