VICAIR® VECTOR 02

THE MOST STABLE CUSHION FOR WHEELCHAIR USERS









PRESSURE ULCER TREATMENT

Nine adjustable compartments for optimal stability

Expanding the user's reach

Extra positioning of the thighs

The Vicair Vector O2 cushion offers enhanced stability to the wheelchair user. This results in improved functionality of the upper body and expands the reach of the user's arms. The cushion's unique compartment layout reduces pressure in the vulnerable area around the sitting bones and tailbone. The Vector O2 provides excellent leg guidance, optimal skin protection and extra pelvic stability.

The cushion is 100% machine washable and breathable, for optimal hygiene and perfect temperature and moisture regulation.

| PRODUCT PERFORMANCE | |
|----------------------------------|-------------|
| Skin protection / Anti decubitus | *** |
| Pressure redistribution | *** |
| Shear force reduction | |
| Microclimate control | *** |
| Impact & vibration damping | AAAA |
| | |

| Positioning | *** |
|---------------------------|-------------|
| Lateral stability | AAAA |
| Frontal stability | AAAA |
| Sliding reduction | AAAA |
| Allows manual adjustments | *** |
| Allows easy transfers | |

| Maintenance & Safety | *** |
|----------------------------|-------------|
| Durability | *** |
| Reliability | *** |
| Functionality when misused | |
| Low maintenance level | AAAA |
| Level of hygiene | *** |
| | |

| Functionality & Application | *** |
|-----------------------------|-------------|
| Comfort | AAAA |
| Low cushion weight | AAAA |
| Available size range | AAAA |
| Maximum user weight | AAAA |
| Low seat to floor height | AAAA |
| Auto adjustments | |

USED WITH

Spinal cord injuries, neuromuscular diseases, hypermobility syndrome, Ehlers Danlos Syndrome (EDS)

| TECHNICAL SPECIFICATIONS | | | | | | | |
|--------------------------|-------------|--|--|--|--|--|--|
| Height | 6 and 10 cm | | | | | | |
| Average weight | 725 grams | | | | | | |
| Max. user weight | 250 kg | | | | | | |
| Machine washable | 60°C max. | | | | | | |

| CL | JSH | 101 | I SI | ZES | 5 | | | | | | | |
|------|----------|-----|------|-----|----|----|-------------------------------|----|----|----|----|-------|
| | Standard | | | | | | (Sizes >53 cm height 10 only) | | | | | |
| | On | req | ues | st | | | | | | | | |
| | | | | | | | | | | | | width |
| Inch | | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 |
| | СМ | 35 | 38 | 40 | 43 | 45 | 48 | 50 | 53 | 55 | 58 | 60 |
| 14 | 35 | | | | | | | | | | | |
| 15 | 38 | | | | | | | | | | | |
| 16 | 40 | | | | | | | | | | | |
| 17 | 43 | | | | | | | | | | | |
| 18 | 45 | | | | | | | | | | | |
| 19 | 48 | | | | | | | | | | | |
| 20 | 50 | | | | | | | | | | | |
| 21 | 53 | | | | | | | | | | | |
| 22 | 55 | | | | | | | | | | | |
| 23 | 58 | | | | | | | | | | | |
| 24 | 60 | | | | | | | | | | | |